



# Safety Funnels

for safe disposal of waste fluids



## Safety funnels



SCAT Safety Funnels open solely when disposing fluids. The ball valve closes automatically and keeps your container reliably closed. The screw cap can rotate freely for proper sealing.



Please find all safety funnels available on next page



Example of use: Safety Funnel with electroconductive PE-HD container



The stainless steel sieve offers:  
- additional splash guard  
- filter for your magnetic stirrers



Please find more safety funnels on pages 17 and 25.



Safety Funnel made of non-conductive PE-HD (without sieve)



## Safety Funnels (PE-HD, electroconductive)



Available with and without sieve. Delivered with antistatic cable.

Part No.	Cap	Material	Part No. with sieve
107672	GL 45	PE-HD (electroconductive), black	108972
107674	S 51	PE-HD (electroconductive), black	108974
107675	S 55	PE-HD (electroconductive), black	108975
107671	S 60 / 61	PE-HD (electroconductive), black	108971
107676	S 65	PE-HD (electroconductive), black	108976
107678	S 70 / 71	PE-HD (electroconductive), black	108978
107677	B 83	PE-HD (electroconductive), black	108977
107673	S 90	PE-HD (electroconductive), black	108973

## Safety Funnels (PE-HD)



Available with and without sieve.

Part No.	Cap	Material	Part No. with sieve
107972	GL 45	PE-HD, white	108672
107974	S 51	PE-HD, white	108674
107975	S 55	PE-HD, white	108675
107971	S 60 / 61	PE-HD, white	108671
107976	S 65	PE-HD, white	108676
107978	S 70 / 71	PE-HD, white	108678
107979	B 83	PE-HD, white	108679
107973	S 90	PE-HD, white	108673

## Stainless Steel Sieve (singular) for Safety Funnels



Part No.	Description
107715	Stainless steel sieve, singular, for safety funnels. Housing: PE-HD, electroconductive, black. Sieve: Stainless steel.



THE SAFER  
SOLUTION!

