



Model	ZHWY-3122B	ZHWY-3222
Control	P.I.D microprocessor	P.I.D microprocessor
Shaking Speed	30-300 rpm (± 1 rpm)	30-300 rpm (± 1 rpm)
Diameter of Orbit	35mm	35mm
Timer	0-500 hours	0-500 hours
Capacity	100mlx136 or 250mlx92 500mlx66 or 1000mlx36	100mlx168 or 250mlx100 500mlx76 or 1000mlx44
Platform Size	800x600 mm	920x640 mm
Platform included	2	2
Exterior Dimension	900x640x780 mm	1100x760x780 mm
Voltage	220V 50/60Hz	220V 50/60Hz
Weight	340Kg	360Kg
Approval	CE, ISO	CE, ISO
Warranty	1 year	1 year