

ZHICHENG[®]


ZHWY-2102C/1102C

Two Layers Incubator Shakers



For more information, please contact our distributor in South Africa: **MICRO FILTRATION CC**
P.O. BOX 334, KLOOF 3640 REP. OF SOUTH AFRICA Contact: Mr. Derick Peters
Tel: +27-31-7621394 FAX +27-31-7633862 Email: microfil@global.co.za

Control Panel



Specifications:

Commodity Name/Model	ZHWY-2102C	ZHWY-1102C
Performance	Heated & Refrigerated	Heated only
Temperature Range	4-60°C	30-60°C
Temperature Accuracy	±0.1°C	±0.1°C
Temperature Uniformity	±1°C	±1°C
Diameter of Orbit	26mm	26mm
Shaking Speed	30-300 rpm (±1rpm)	30-300 rpm (±1rpm)
Timer	1 minute to 500 hours	1 minute to 500 hours
Standard Capacity	250mlx30	250mlx30
Max.Capacity	100mlx48 250mlx30 500mlx24 50mlx48	100mlx48 250mlx30 500mlx24 50mlx48
Platform Size	500x350 mm	500x350 mm
Platform included	2	2
Exterior Dimension	700x590x1230 mm	700x590x1230 mm
Motor Power	550W	436W
Electrical Requirement	220V 50/60Hz	220V 50/60Hz
Net Weight	165Kg	150Kg
Gross Weight	190Kg	175Kg
Manufacturer	ZHICHENG (CHINA)	ZHICHENG (CHINA)
Warranty	1 year	1 year
Approval	CE, ISO	CE, ISO

ZHICHENG[®] Large Capacity Two Layers Incubator Shakers



Model	ZHWY-2102	ZHWY-1102	ZHWY-2112B	ZHWY-1112B
Performance	Heated & Refrigerated	Heated only	Heated & Refrigerated	Heated only
Temperature Range	4-60°C	30-60°C	4-60°C	30-60°C
Temperature Accuracy	±0.1°C	±0.1°C	±0.1°C	±0.1°C
Temperature Uniformity	±1°C	±1°C	±1°C	±1°C
Diameter of Orbit	26mm	26mm	35mm	35mm
Shaking Speed	30-300 rpm(±1 rpm)	30-300 rpm(±1 rpm)	30-300 rpm(±1 rpm)	30-300 rpm(±1 rpm)
Timer	0-500 hours	0-500 hours	0-500 hours	0-500 hours
Capacity	100mlx90 250mlx56 500mlx48 1000mlx30	100mlx90 250mlx56 500mlx48 1000mlx30	100mlx154 250mlx90 500mlx70 1000mlx36	100mlx154 250mlx90 500mlx70 1000mlx36
Platform Size	734x458 mm	734x458 mm	970x560mm	970x560mm
Platform included	2	2	2	2
Exterior Dimension	950x710x1250 mm	950x710x1250 mm	1200x820x1600 mm	1200x820x1600 mm
Motor Power	560W	446W	961W	446W
Voltage	220V 50/60Hz	220V 50/60Hz	220V 50/60Hz	220V 50/60Hz
Weight	272Kg	265Kg	450Kg	446Kg
Warranty	1 year	1 year	1 year	1 year
Approval	CE	CE	CE	CE

For more information, please contact our distributor in South Africa: **MICRO FILTRATION CC**
P.O. BOX 334, KLOOF 3640 REP. OF SOUTH AFRICA Contact: Mr. Derick Peters
Tel: +27-31-7621394 FAX +27-31-7633862 Email: microfil@global.co.za